

Active[®]

by ACCESS[™]



Active Cane holder installation

Active




1.
Caneholder kit
- Cane cup
 - Fastening brackets
 - Cane clamp
 - 4 screws

Tool needed

- Screwdriver (PH2)



2.
Fasten the cane cup and the cane clamp in the positions shown on the picture. use the brackets and the enclosed screws.



Kuhn und Bieri AG

Rehaprodukte und Dienstleistungen für Heime, Spitäler und Spitex
Produits réha et services pour homes, hôpitaux et organismes de soins à domicile

Sägestrasse 75, CH-3098 Köniz
Telefon 0848 102040, Fax 031 970 01 71, info@kuhnbieri.ch

www.kuhnbieri.ch